

Visium Total: spatial total RNA sequencing via in situ polyA-tailing

Last Edits Made: DWM, 6/29/2021

Materials/Reagents

- Yeast PAP
- 5X Yeast PAP reaction buffer
- 20U/ul SUPERase in RNase Inhibitor

Buffers

Yeast PAP Reaction Mix:

	Final Conc.	Volume to add			
		1 RXN	4.4 RXNs	8.8 RXNs	
5X Yeast PAP reaction buffer	1X	15ul	66ul	132ul	
Yeast PAP	0.25M	3ul	13.2ul	26.4ul	
25mM ATP	1%	1.5ul	6.6ul	13.2ul	
SUPERase	0.8U/ul	3ul	13.2ul	26.4ul	
H ₂ O	--	52.5ul	231ul	462ul	
Total Volume	--	75ul	330ul	660ul	

Yeast Reaction Buffer:

	Final Conc.	Volume to add			
		1 RXN	4.4 RXNs	8.8 RXNs	
5X Yeast PAP reaction buffer	1X	50ul			
SUPERase	0.2U/ul	2.5ul			
H ₂ O	--	197.5ul			
Total Volume	--	250ul			

Protocol

1. Section onto Visium slide according to protocols from 10x Genomics.
2. Fix, stain, and image tissue according to Visium protocol ("Step 1: Tissue Staining & Imaging").
3. Prepare appropriate volumes of **Yeast PAP Reaction Mix** and **Yeast Reaction Buffer**.
4. Set the thermocycler to 37°C (Step 2.1.a)
5. Place the slide into the Slide Cassette (See Visium protocol for "Tips and Best Practices").
6. Add **75ul** of **Yeast Reaction Buffer** to each well to equilibrate tissue. Incubate at room temperature for 30sec.
7. Remove and discard buffer from each well.
8. Add **75ul** of **Yeast PAP Enzyme Mix**.
9. Place slide onto the Thermocycler Adaptor and incubate at 37°C for 25min.
10. Remove and discard enzyme mix from each well. Proceed directly to Step 2.1.c to permeabilize tissue, and follow the standard Visium protocol to completion.